BBrain Basic Clock

NEDERLANDS **3**

ENGLISH **13**

DEUTSCH **22**

FRANÇAIS **32**

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10

11

Thank you for choosing the **BBrain clock**! This calendar clock supports seniors, individuals with visual impairments, and those with cognitive challenges such as dementia or acquired brain injury, in maintaining orientation in time and daily planning and structure. This manual explains how the **BBrain clock** works and how to set it up according to your personal preferences.


1. About this manual

- Before using the **BBrain Clock**, carefully read this manual, including the safety information, ensuring proper and safe handling of the device.
- Images and screenshots in this manual may differ from the actual product.
- Available functions and services depend on the device's hardware and software.
- BBrain is not liable for compatibility issues caused by user modifications or bad usage.
- All audio, images, and backgrounds on this device are licensed only for personal use. Using them for commercial purposes is a violation of copyright laws. BBrain is not responsible for any copyright infringements by users.

2. Package contents

Before proceeding with the installation, ensure all necessary parts are included. If any items are missing, please contact us at: info@bbrain.eu.

- 1 BBrain Calendar Clock
- 1 BBrain Stand
- 1 Charger
with a 2-meter cable
- 1 Remote Controller
(battery included)
- 1 User Manual


3. Startup

The **BBrain Clock** will be operational after following the next initial steps:

- Plug the energized charger cable into the device, and it will turn on automatically.
- The startup screen will be displayed for approximately 5 seconds (1), followed by the main time interface (2).

1

2

4. Setting up the calendar clock

The following steps explain how to interact with the **Bbrain Clock**:

- Press the **"MENU"** button on the back of the clock or the **"EXIT/SETUP"** button on the remote controller to access the settings.
- Use the  and  -buttons to navigate through the settings of the menu. The selected feature will change to a yellow color.
- Directly use the  and  -buttons to modify the state of each feature.
- If you want to edit the value of a feature, press the **"OK"** button on the back of the clock or the **"OK/ENTER"** button on the remote controller to start editing, then use the  and  -buttons to select the desired values and finish the editing process by pressing the **"MENU"** button on the back of the clock or the **"EXIT/SETUP"** button on the remote controller to save the selected value. When multiple options are available, the selected values will be displayed in blue color.
- Press the **"MENU"** button on the back of the clock or the **"EXIT/SETUP"** button on the remote controller to save changes and return to the main clock screen.

5. Features description

This section lists all the necessary features to complete the initial setup of the **Bbrain Clock**:

LANGUAGE	< ENGLISH >
----------	-------------

- Select your preferred language: more than 10 options available.
Directly use the  and  -buttons to modify the desired language.

TIME	< 00:00 >
------	-----------

- Set the current local time.
Press **"OK"** to set the time.
Adjust hours and minutes using  and  -buttons, with the selected field underlined.
Confirm the time by pressing **"OK"**.

Note: Select the 24-hour mode before setting the time to ensure the correct AM/PM display.

DATE	< DD-MM-JJJJ >
------	----------------

- Set the actual date.
Press **"OK"** to set the date. Use the  and  -buttons to select and adjust the day, month, or year (e.g., DD/MM/YYYY). Confirm your settings by pressing **"OK"**.

TIME MODE	< 24 >
-----------	--------

- Choose between 12-hour (AM/PM) or 24-hour format.
Directly use the  and  -buttons to modify the time mode.

DATE MODE	< Day - Month - Year >
-----------	------------------------

- Choose a preferred date format:
Day/Month/Year
Month/Day/Year
Year/Month/Day
Directly use the  and  -buttons to modify the date mode.

DAY PERIODS	< 06 >
--------------------	---------------------

- For all languages, the day will include four distinct parts, allowing users to customize the start time for each segment.

THE MORNING STARTS AT

THE AFTERNOON STARTS AT

THE EVENING STARTS AT

THE NIGHT STARTS AT

- For Flemish and German languages, the day includes 2 more parts:

THE PRE-AFTERNOON STARTS AT

THE POST-AFTERNOON STARTS AT

Directly use the ◀ and ▶ -buttons to modify the hour when each part of the day starts.

SELECT HOUR FOR TONE	< 1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16 17 18 19 20 21 22 23 00 >
-----------------------------	---

- You can set a tone as an additional reminder on specific hours.
Use the ◀ and ▶ -buttons to select an hour (highlighted with a white underline).
Press **"OK"** to activate the tone. Programmed hours will appear in blue.
To delete or modify a reminder, select it and press **"OK"**.

BRIGHTNESS DAY	< 15 >
-----------------------	---------------------

- Set the daytime brightness.
Directly use the ◀ and ▶ -buttons to modify the desired brightness.

BRIGHTNESS NIGHT	< 01 >
-------------------------	---------------------

- Set the nighttime brightness.
Directly use the ◀ and ▶ -buttons to modify the desired brightness.

NIGHTMODE	< OFF >
------------------	----------------------

- Set the hour after which the screen color changes to black with yellow text, allowing a restful night without distractions.
Directly use the ◀ and ▶ -buttons to modify the desired hour.


VOLUME	< 04 >
--------	--------

- Choose a preferred volume to reproduce the tone and alarms:
Directly use the  and  -buttons to modify the desired volume.

COLOR MODE	< BLACK - WHITE TEXT >
------------	------------------------

- Choose one of three screen color modes as shown (3):

Black/White Text: Standard view.


Black/Yellow Text: Suitable for the visually impaired.

Black/Red Bar: Intended to capture the user's attention.

Directly use the  and  -buttons to select the desired color mode.

SCREEN INTERVAL MODE	< OFF >
----------------------	---------

- The **BBrain Clock** has 2 different screen interfaces to show the date and time information. One complete (4) and one basic (5):

You can manually switch the interface by pressing and holding the “OK” button on the back of the clock for 3 seconds.

The screen interval mode allows users to either maintain a consistent interface or cycle between the two available options at intervals of 30, 60, or 90 seconds.

- Directly use the  and  -buttons to select the desired interval mode.

ALARM SETTINGS	REPEAT < 1 >	DURATION < 30 >	>
----------------	--------------	-----------------	---

- The **BBrain Clock** allows six different customizable alarms that can be programmed for specific dates, providing ideal reminders for daily routines and important holidays.
Press **"OK"** to set the alarms.
Directly use the ◀ and ▶ -buttons to select the parameter to modify.
Use the ◀ and ▶ -buttons to modify the Repeat and Duration values to apply for all the alarms:
 - o Configure the repetition intervals as 0, 1, 2, or 3, using a 5-minute cycle if the user does not press any button to acknowledge the alarm while it is active.
 - o Configure the tone duration to 30, 60, or 90 seconds. Press the **"OK"** button while the > icon is highlighted to enter the Advanced Alarm Menu.
 - o Use the ▲ and ▼ -buttons to scroll through the following alarm options and press **"OK"** to confirm your selection:

WAKE UP ALARM

GO TO BED ALARM

TAKE A NAP ALARM


USER'S BIRTHDAY ALARM

CHRISTMAS DAY ALARM

NEW YEAR'S DAY ALARM

- o Use the ◀ and ▶ -buttons to navigate to the desired day(s) of the week. Active days are highlighted with a white underline.
- o Press **"OK"** to activate the alarm for the selected day. Active days will appear in blue.
- o Adjust the hours and minutes using the ▲ and ▼ -buttons. The selected field will be underlined.
- o Press **"OK"** to confirm the time.
- o Press the **"MENU"** button to exit the Advanced Alarm Menu.

Below are some examples illustrating the appearance of the alarms:

MEDICATION SETTINGS	REPEAT < 1 >	DURATION < 30 >	>
---------------------	--------------	-----------------	---

- The **BBrain Clock** features three customizable medication alarms. Press **"OK"** to set the medicine alarms.

Directly use the and -buttons to select the parameter to modify.

Use the and -buttons to modify the Repeat and Duration values to apply for all the medicine alarms:

- o Configure the repetition intervals as 0, 1, 2, or 3, using a 5-minute cycle if the user does not press any button to acknowledge the alarm while it is active.
- o Configure the tone duration to 30, 60, or 90 seconds.
- o Press **OK** button while the **>** icon is highlighted to enter the Advanced Medication Alarms Menu.
- o Use the and -buttons to scroll through the following alarm options and press **"OK"** to confirm your selection:

MORNING MEDICATION

AFTERNOON MEDICATION

EVENING MEDICATION

- o Use the and -buttons to navigate to the desired day(s) of the week. Active days are highlighted with a white underline.
- o Press **"OK"** to activate the medication alarm for the selected day. Active days will appear in blue.
- o Repeat this step for all desired days.
- o Adjust the hours and minutes using the and -buttons. The selected field will be underlined.
- o Press **"OK"** to confirm the time.
- o Press the **"MENU"** button to exit the Advanced Medication Alarm Menu.

Below are some examples illustrating the appearance of the medication alarms:


RESET SETTINGS	< No >
----------------	--------

- Allows the user to restore the **BBrain Clock** to factory default settings. Only the set “**DATE**” and “**TIME**” set will remain without changes, the rest of the settings, including the alarms, will be lost. Use the ◀ and ▶ -buttons to select the desired option. A popup message will appear asking for confirmation to Reset the settings (6).

6

Use the ◀ and ▶ -buttons to select “**YES**” or “**NO**”. The selected field will appear in yellow. Confirm the desired option by pressing the “**OK**” button.

FIRMWARE VERSION	SW4. XXXXXX
------------------	-------------

- The current software version is displayed for internal control purpose.

6. Safety instructions and maintenance

It is crucial to adopt the following key recommendations to ensure optimal performance and significantly extend the product’s lifespan.

6.1 Indoor Use Only:

The **BBrain Clock** is designed exclusively for indoor use. Do not attempt to operate the device or any of its accessories outdoors under any circumstances.

6.2 Power source:

The **BBrain Clock** must be connected to a standard household power source. Connecting it to unregulated or unauthorized power sources for domestic use may result in damage to the device

- 6.3 Use the **BBrain Clock** exclusively with the provided power adapter.

- 6.4 **Placement Recommendation:**
Position the **BBrain Clock** on a piece of furniture that provides a clear, direct view for the user throughout most of the day. Make sure it is securely placed to prevent falls.
- 6.5 **Avoiding Screen Damage:**
Avoid scratching the screen with any sharp or pointed objects.
- 6.6 **Operating Temperature and Humidity:**
Operate the **BBrain Clock** within a temperature range of 10°C to 40°C (50°F to 104°F) and a humidity range of 20% to 80% (non-condensing) for optimal performance.
- 6.7 **Cleaning Instructions:**
Use a soft, dry cloth to clean the surface and screen of the **BBrain Clock**. Avoid using abrasive materials, liquids, or chemical cleaners.
- 6.8 **Remote Controller Battery Management:**
If the remote controller stops working, it is recommended to replace the CR2025 3V Lithium Manganese Dioxide coin battery as shown in the image of its back.

7. BBrain clock pictures

BBrain Clock: reference images of its front (7) and back view (8)

7

8

Remote Controller: reference images of its front (9), rear (10) and compatible battery (11)

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